

Product datasheet for **TA180014M**

Mouse monoclonal DsRed-Express2 antibody, clone OT11C4


Product data:

Product Type:	Tag Antibodies
Product Name:	Mouse monoclonal DsRed-Express2 antibody, clone OT11C4
Clone Name:	OT11C4 (formerly 1C4)
Applications:	WB
Recommended Dilution:	WB: 1:2000
Reactivity:	DsRed-Express2
Host:	Mouse
Clonality:	Monoclonal
Immunogen:	This antibody is raised against a recombinant protein of DsRed-Express2 expressed in HEK293 cells.
Isotype:	IgG2b
Formulation:	PBS (pH 7.3) containing 1%BSA, 50% glycerol and 0.02% sodium azide
Concentration:	3.64mg/ml
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	Stable for 1 year at -20°C from date of shipment.
Predicted Protein Size:	25.7 kDa



[View online »](#)

Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY DsRed-Express2 (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-DsRed-Express2 ([TA180014], 1:2000).