


Product datasheet for TA180062S

Mouse Monoclonal mOrange Antibody, clone OTI2B11

Product data:

Product Type:	Tag Antibodies
Product Name:	Mouse Monoclonal mOrange Antibody, clone OTI2B11
Clone Name:	OTI2B11
Applications:	IHC, WB
Recommended Dilution:	WB 1:2000, IHC 1:100
Reactivity:	mOrange
Host:	Mouse
Clonality:	Monoclonal
Immunogen:	A mOrange tagged recombinant protein expressed in E.coli.
Isotype:	IgG2a
Formulation:	PBS (pH7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide
Concentration:	1 mg/ml
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	Stable for 1 year at -20°C from date of shipment.


Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY mOrange (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-mOrange.



[View online »](#)



Immunohistochemical staining of Frozen section Brain tissue from hGFAP-Cre, Rosa26 TdTomato mice using anti-mOrange mouse monoclonal antibody (1:100). Secondary antibody: Alexa 488 goat anti-mouse.