

Product datasheet for **TA180072S**


Mouse monoclonal mPlum antibody, clone OT12G4

Product data:

Product Type:	Tag Antibodies
Product Name:	Mouse monoclonal mPlum antibody, clone OT12G4
Clone Name:	OT12G4 (formerly 2G4)
Applications:	WB
Recommended Dilution:	WB: 1:2000
Reactivity:	mPlum
Host:	Mouse
Clonality:	Monoclonal
Immunogen:	This antibody is raised against a recombinant protein of mPlum expressed in HEK293 cells.
Isotype:	IgG1
Formulation:	PBS (pH7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide
Concentration:	2.67mg/ml
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	Stable for 1 year at -20°C from date of shipment.
Predicted Protein Size:	25.6 kDa



[View online »](#)

Product images:

HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY mPlum (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-mPlum ([TA180072], 1:2000).