

Product datasheet for **TA180083S**


Mouse monoclonal DsRed-Monomer antibody, clone OTI4C7

Product data:

Product Type:	Tag Antibodies
Product Name:	Mouse monoclonal DsRed-Monomer antibody, clone OTI4C7
Clone Name:	OTI4C7 (formerly 4C7)
Applications:	WB
Recommended Dilution:	WB: 1:2000
Reactivity:	DsRed-Monomer
Host:	Mouse
Clonality:	Monoclonal
Immunogen:	This antibody is raised against a recombinant protein of DsRed-Monomer expressed in HEK293 cells.
Isotype:	IgG1
Formulation:	PBS (pH7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide
Concentration:	1 mg/ml
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	Stable for 1 year at -20°C from date of shipment.
Predicted Protein Size:	25.4 kDa



[View online »](#)

Product images:

HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY DsRed-Monomer (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-DsRed-Monomer ([TA180083], 1:2000).