

## Product datasheet for **TA180092S**


### Mouse monoclonal mBFP antibody, clone OT12E7

#### Product data:

Product Type:	Tag Antibodies
Product Name:	Mouse monoclonal mBFP antibody, clone OT12E7
Clone Name:	OT12E7 (formerly 2E7)
Applications:	WB
Recommended Dilution:	WB: 1:2000
Reactivity:	mBFP
Host:	Mouse
Clonality:	Monoclonal
Immunogen:	A mBFP tagged recombinant protein expressed in HEK293T cell.
Isotype:	IgG1
Formulation:	PBS (pH7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide
Concentration:	1 mg/ml
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	Stable for 1 year at -20°C from date of shipment.
Predicted Protein Size:	27.1 kDa
Note:	This antibody also recognizes mRFP protein.



[View online »](#)

**Product images:**

HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY mRFP (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-mRFP ([TA180092], 1:2000).