

Product datasheet for **TA180128**

HA mouse monoclonal antibody, clone CB051


Product data:

Product Type:	Tag Antibodies
Product Name:	HA mouse monoclonal antibody, clone CB051
Clone Name:	CB051
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	HA tag
Host:	Mouse
Clonality:	Monoclonal
Immunogen:	Monoclonal antibody is produced by immunizing mice with a synthetic peptide from influenza hemagglutinin coupled to KLH
Isotype:	IgG1
Formulation:	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	Stable for 1 year at -20°C from date of shipment.
Synonyms:	ANTI-HA; Anti-HA Antibody; Anti HA; Anti HA Antibody; ANTI INFLUENZA




[View online »](#)

Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or HA tagged LGALS3 cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-HA.



Western Blot comparison of anti-HA Tag antibodies from OriGene (TA180128), BioLegend (mms-101p), and Abcam (ab18181). 0.04ul (Lane 1, 3, 5) and 0.2ul (Lane 2, 4, 6) of HA-PARD6B lysate (recombinant proteins with HA tag at the N-terminus), blotted with 1:1000 dilution of TA180128 (lane 1-2), mms-101p (Lane 3-4), and ab18181 (Lane 5-6).