

Product datasheet for **TA180144S**

DDK (FLAG) mouse monoclonal antibody, clone OTI11C3


Product data:

Product Type:	Tag Antibodies
Product Name:	DDK (FLAG) mouse monoclonal antibody, clone OTI11C3
Clone Name:	OTI11C3
Applications:	ICC, IF, IHC, WB
Recommended Dilution:	WB 1:4000
Reactivity:	DDK (FLAG tag): DYKDDDDK
Host:	Mouse
Clonality:	Monoclonal
Immunogen:	Anti-DDK monoclonal antibody is produced by immunizing mice with a synthetic peptide coupled to KLH
Isotype:	IgG1
Formulation:	PBS (pH7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide
Concentration:	1 mg/ml
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	Stable for 1 year at -20°C from date of shipment.
Note:	This anti-DDK (FLAG) antibody recognizes both N-terminus and C-terminus DDK tags. It is recommended for Western blot only.




[View online »](#)

Product images:



HEK293T cells were transfected with the pCMV6-ENTRY YPEL3 ([RC208686]; Lane 1), CYSLTR1 ([RC207763]; Lane 2), SEMA3D ([RC216032]; Lane 3), Dmp1 ([RC222208]; Lane 4) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-DDK (1:4000).



Immunofluorescent staining of 293T cells transfected by human CTAG1B plasmid([RC213318]) using anti-DDK antibody ([TA180144]/green, right) compared to an Isotype control(Left).(1:400)