

Product datasheet for TA180148S

Mouse monoclonal tdTomato antibody, Clone OTI4B3

Product data:

Product Type: Tag Antibodies

Product Name: Mouse monoclonal tdTomato antibody, Clone OTI4B3

Clone Name: OTI4B3

Applications: WB

Reactivity: WB: 1:500
Host: WB: 1:500

WB: 1:500

Mouse

Clonality: Monoclonal

Immunogen: A tdTomato tagged recombinant protein expressed in E.coli.

Isotype: lgG⁻

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Conjugation: Unconjugated Storage: Store at -20°C.

Stability: Stable for 1 year at -20°C from date of shipment.

Predicted Protein Size: 54.2 kDa


OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com

EU: info-de@origene.com CN: techsupport@origene.cn



Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY tdTomato (JB1016_D07, Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-tdTomato (1:500).