

Product datasheet for **TA183001**

Mouse monoclonal ZsYellow1 antibody, HRP conjugated

Product data:

Product Type:	Tag Antibodies
Product Name:	Mouse monoclonal ZsYellow1 antibody, HRP conjugated
Clone Name:	OT13C4 (formerly 3C4)
Applications:	WB
Recommended Dilution:	WB: 1:1000-5000
Reactivity:	ZsYellow1
Host:	Mouse
Clonality:	Monoclonal
Immunogen:	This antibody is raised against a recombinant protein and recognizes ZsYellow1.
Isotype:	IgG1
Formulation:	PBS (pH7.4) containing 0.1% BSA and 50% glycerol
Concentration:	0.5mg/ml
Conjugation:	HRP
Storage:	Store at -20°C. For maximum recovery of product, centrifuge the original vial after thawing and prior to opening the cap. Aliquot to avoid repeated freezing and thawing.
Stability:	Stable for 1 year at -20°C from date of shipment.
Predicted Protein Size:	26 kDa



[View online »](#)