

## Product datasheet for **TA183002**

### Mouse monoclonal DsRed2 antibody, HRP conjugated

#### Product data:

Product Type:	Tag Antibodies
Product Name:	Mouse monoclonal DsRed2 antibody, HRP conjugated
Clone Name:	OT11B11 (formerly 1B11)
Applications:	WB
Recommended Dilution:	WB: 1:1000-5000
Reactivity:	DsRed2
Host:	Mouse
Clonality:	Monoclonal
Immunogen:	This antibody is raised against a recombinant protein and recognizes the DsRed2 protein.
Isotype:	IgG2b
Formulation:	PBS (pH7.4) containing 0.1% BSA and 50% glycerol
Concentration:	0.5mg/ml
Conjugation:	HRP
Storage:	Store at -20°C. For maximum recovery of product, centrifuge the original vial after thawing and prior to opening the cap. Aliquot to avoid repeated freezing and thawing.
Stability:	Stable for 1 year at -20°C from date of shipment.
Predicted Protein Size:	26 kDa



[View online »](#)