

Product datasheet for TA183004

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Mouse monoclonal tdTomato antibody, HRP conjugated

Product data:

Product Type: Tag Antibodies

Product Name: Mouse monoclonal tdTomato antibody, HRP conjugated

Clone Name: OTI2H2 (formerly 2H2)

Applications: WB

Recommended Dilution: WB: 1:1000-5000

Reactivity: tdTomato
Host: Mouse

Clonality: Monoclonal

Immunogen: This antibody is raised against a recombinant protein and recognizes tdTomato.

Isotype: IgG2b

Formulation: PBS (pH7.4) containing 0.1% BSA and 50% glycerol

Concentration: 0.5mg/ml

Conjugation: HRP

Storage: Store at -20°C. For maximum recovery of product, centrifuge the original vial after thawing


and prior to opening the cap. Aliquot to avoid repeated freezing and thawing.

Stability: Stable for 1 year at -20°C from date of shipment.

Predicted Protein Size: 54.2 kDa



Product images:



HEK293T cells were transfected with pCMV6-ENTRY-tdTomato cDNA for 48 hrs and lysed. Cell lysate (10 ug) was separated by SDS-PAGE and immunoblotted with HRP conjugated antitdTomato monoclonal antibody (TA183004, 1:1000)