

Product datasheet for **TA183006**

Mouse monoclonal DsRed-Express2 antibody, HRP conjugated

Product data:

| | |
|-------------------------|---|
| Product Type: | Tag Antibodies |
| Product Name: | Mouse monoclonal DsRed-Express2 antibody, HRP conjugated |
| Clone Name: | OT11C4 (formerly 1C4) |
| Applications: | WB |
| Recommended Dilution: | WB: 1:1000-5000 |
| Reactivity: | DsRed-Express2 |
| Host: | Mouse |
| Clonality: | Monoclonal |
| Immunogen: | This antibody is raised against a recombinant protein of DsRed-Express2 expressed in HEK293 cells. |
| Isotype: | IgG2b |
| Formulation: | PBS (pH7.4) containing 0.1% BSA and 50% glycerol |
| Concentration: | 0.5mg/ml |
| Conjugation: | HRP |
| Storage: | Store at -20°C. For maximum recovery of product, centrifuge the original vial after thawing and prior to opening the cap. Aliquot to avoid repeated freezing and thawing. |
| Stability: | Stable for 1 year at -20°C from date of shipment. |
| Predicted Protein Size: | 25.7 kDa |



[View online »](#)