

Product datasheet for TA183012

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Mouse monoclonal mOrange/mOrange2 antibody, HRP conjugated

Product data:

Product Type: Tag Antibodies

Product Name: Mouse monoclonal mOrange/mOrange2 antibody, HRP conjugated

Clone Name: OTI4E10 (formerly 4E10)

Applications: WB

Recommended Dilution: WB: 1:1000-5000

Reactivity: mOrange, mOrange2

Host: Mouse

Clonality: Monoclonal

Immunogen: This antibody is raised against a recombinant protein of mOrange expressed in HEK293 cells.

It recognizes both mOrange and mOrange2 proteins.

Isotype: IgG2b

Formulation: PBS (pH7.4) containing 0.1% BSA and 50% glycerol

Concentration: 0.5mg/ml

Conjugation: HRP

Storage: Store at -20°C. For maximum recovery of product, centrifuge the original vial after thawing

and prior to opening the cap. Aliquot to avoid repeated freezing and thawing.

Stability: Stable for 1 year at -20°C from date of shipment.

Predicted Protein Size: 26.7 kDa