

Product datasheet for **TA183015**

Mouse monoclonal mGFP/mCFP/mYFP/KillerRed antibody, HRP conjugated

Product data:

Product Type:	Tag Antibodies
Product Name:	Mouse monoclonal mGFP/mCFP/mYFP/KillerRed antibody, HRP conjugated
Clone Name:	OT12F6 (formerly 2F6)
Applications:	WB
Recommended Dilution:	WB: 1:1000-5000
Reactivity:	mGFP, mCFP, mYFP, KillerRed
Host:	Mouse
Clonality:	Monoclonal
Immunogen:	This antibody is raised against a recombinant protein of mGFP expressed in HEK293 cells. It recognizes mGFP, mCFP, mYFP and KillerRed proteins.
Isotype:	IgG1
Formulation:	PBS (pH7.4) containing 0.1% BSA and 50% glycerol
Concentration:	0.5mg/ml
Conjugation:	HRP
Storage:	Store at -20°C. For maximum recovery of product, centrifuge the original vial after thawing and prior to opening the cap. Aliquot to avoid repeated freezing and thawing.
Stability:	Stable for 1 year at -20°C from date of shipment.
Predicted Protein Size:	27.3 kDa



[View online »](#)