

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## Product datasheet for TA183016

## Mouse monoclonal DsRed-Monomer antibody, HRP conjugated

## **Product data:**

Product Type:	Tag Antibodies
Product Name:	Mouse monoclonal DsRed-Monomer antibody, HRP conjugated
Clone Name:	OTI4C7 (formerly 4C7)
Applications:	WB
Recommended Dilution:	WB: 1:1000-5000
Reactivity:	DsRed-Monomer
Host:	Mouse
Clonality:	Monoclonal
Immunogen:	This antibody is raised against a recombinant protein of DsRed-Monomer expressed in HEK293 cells.
lsotype:	lgG1
Formulation:	PBS (pH7.4) containing 0.1% BSA and 50% glycerol
Concentration:	0.5mg/ml
Conjugation:	HRP
Storage:	Store at -20°C. For maximum recovery of product, centrifuge the original vial after thawing and prior to opening the cap. Aliquot to avoid repeated freezing and thawing.
Stability:	Stable for 1 year at -20°C from date of shipment.
Predicted Protein Size:	25.4 kDa



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US