

Product datasheet for **TA183022**

Mouse monoclonal AsRed2 antibody, HRP conjugated

Product data:

Product Type:	Tag Antibodies
Product Name:	Mouse monoclonal AsRed2 antibody, HRP conjugated
Clone Name:	OT12E2 (formerly 200)
Applications:	WB
Recommended Dilution:	WB: 1:1000-5000
Reactivity:	AsRed2
Host:	Mouse
Clonality:	Monoclonal
Immunogen:	A AsRed2 tagged recombinant protein expressed in HEK293T cell.
Isotype:	IgG1
Formulation:	PBS (pH7.4) containing 0.1% BSA and 50% glycerol
Concentration:	0.5mg/ml
Conjugation:	HRP
Storage:	Store at -20°C. For maximum recovery of product, centrifuge the original vial after thawing and prior to opening the cap. Aliquot to avoid repeated freezing and thawing.
Stability:	Stable for 1 year at -20°C from date of shipment.
Predicted Protein Size:	25.8 kDa



[View online »](#)