

## **Product datasheet for TA183023**

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## Mouse Monoclonal ZsGreen1 Antibody, Clone OTI2C2, Magnetic Beads

## **Product data:**

**Product Type:** Tag Antibodies

**Product Name:** Mouse Monoclonal ZsGreen1 Antibody, Clone OTI2C2, Magnetic Beads

Clone Name: OTI2C2 (formerly 2C2)

Applications: IP, Protein Purification

Recommended Dilution: IP: Usage Dependent

**Reactivity:** ZsGreen1 Host: Mouse

Clonality: Monoclonal

**Immunogen:** This antibody is raised against a recombinant protein of ZsGreen1 expressed in HEK293 cells.

**Isotype:** IgG1

Formulation: PBS (pH7.4) containing 0.1% BSA, 0.02% NaN3 and 50% glycerol with final concentration

5mg/mL beads

**Conjugation:** Magnetic Beads

Storage: Shipped at -20°C or with ice packs, Upon delivery store at -20°C. Dilute in PBS(pH7.3) if

necessary. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws.

Stability: Stable for 1 year at -20°C from date of shipment.

**Predicted Protein Size:** 26 kDa