

## **Product datasheet for TA183025S**

## **OriGene Technologies, Inc.** 9620 Medical Center Drive, Ste 200

Rockville, MD 20850, US
Phone: +1-888-267-4436
https://www.origene.com
techsupport@origene.com
EU: info-de@origene.com
CN: techsupport@origene.cn

## Mouse Monoclonal DsRed2 Antibody, Magnetic Beads

## **Product data:**

**Product Type:** Tag Antibodies

**Product Name:** Mouse Monoclonal DsRed2 Antibody, Magnetic Beads

Clone Name: OTI1B11 (formerly 1B11)

Applications: IP, Protein Purification

Recommended Dilution: IP: Usage Dependent

Reactivity: DsRed2
Host: Mouse

Clonality: Monoclonal

**Immunogen:** This antibody is raised against a recombinant protein and recognizes the DsRed2 protein.

**Isotype:** lgG2b

Formulation: PBS (pH7.4) containing 0.1% BSA, 0.02% NaN3 and 50% glycerol with final concentration

5mg/mL beads

**Conjugation:** Magnetic Beads

Storage: Shipped at -20°C or with ice packs, Upon delivery store at -20°C. Dilute in PBS(pH7.3) if

necessary. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws.

**Stability:** Stable for 1 year at -20°C from date of shipment.

Predicted Protein Size: 26 kDa