

Product datasheet for TA183027S

Mouse Monoclonal tdTomato Antibody, Magnetic Beads

Product data:

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product Type:	Tag Antibodies
Product Name:	Mouse Monoclonal tdTomato Antibody, Magnetic Beads
Clone Name:	OTI2H2 (formerly 2H2)
Applications:	IP, Protein Purification
Recommended Dilution:	IP: Usage Dependent
Reactivity:	tdTomato
Host:	Mouse
Clonality:	Monoclonal
Immunogen:	This antibody is raised against a recombinant protein and recognizes tdTomato.
lsotype:	lgG2b
Formulation:	PBS (pH7.4) containing 0.1% BSA, 0.02% NaN3 and 50% glycerol with final concentration 5mg/mL beads
Conjugation:	Magnetic Beads
Storage:	Shipped at -20°C or with ice packs, Upon delivery store at -20°C. Dilute in PBS(pH7.3) if necessary. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws.
Stability:	Stable for 1 year at -20°C from date of shipment.
Predicted Protein Size:	54.2 kDa



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US