

## Product datasheet for **TA183037S**

### Mouse Monoclonal mPlum Antibody, Magnetic Beads

#### Product data:

Product Type:	Tag Antibodies
Product Name:	Mouse Monoclonal mPlum Antibody, Magnetic Beads
Clone Name:	OT12G4 (formerly 2G4)
Applications:	IP, Protein Purification
Recommended Dilution:	IP: Usage Dependent
Reactivity:	mPlum
Host:	Mouse
Clonality:	Monoclonal
Immunogen:	This antibody is raised against a recombinant protein of mPlum expressed in HEK293 cells.
Isotype:	IgG1
Formulation:	PBS (pH7.4) containing 0.1% BSA, 0.02% NaN <sub>3</sub> and 50% glycerol with final concentration 5mg/mL beads
Conjugation:	Magnetic Beads
Storage:	Shipped at -20°C or with ice packs, Upon delivery store at -20°C. Dilute in PBS(pH7.3) if necessary. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws.
Stability:	Stable for 1 year at -20°C from date of shipment.
Predicted Protein Size:	25.6 kDa



[View online »](#)