

Product datasheet for TA183038

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Mouse Monoclonal mGFP/mCFP/mYFP/KillerRed Antibody, Magnetic Beads

Product data:

Product Type: Tag Antibodies

Product Name: Mouse Monoclonal mGFP/mCFP/mYFP/KillerRed Antibody, Magnetic Beads

Clone Name: OTI2F6 (formerly 2F6)

Applications: IP, Protein Purification

Recommended Dilution: IP: Usage Dependent

Reactivity: mGFP, mCFP, mYFP, KillerRed

Host: Mouse

Clonality: Monoclonal

Immunogen: This antibody is raised against a recombinant protein of mGFP expressed in HEK293 cells. It

recognizes mGFP, mCFP, mYFP and KillerRed proteins.

Isotype: lgG1

Formulation: PBS (pH7.4) containing 0.1% BSA, 0.02% NaN3 and 50% glycerol with final concentration

5mg/mL beads

Conjugation: Magnetic Beads

Storage: Shipped at -20°C or with ice packs, Upon delivery store at -20°C. Dilute in PBS(pH7.3) if

necessary. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws.

Stability: Stable for 1 year at -20°C from date of shipment.

Predicted Protein Size: 27.3 kDa

