

Product datasheet for **TA183038S**

Mouse Monoclonal mGFP/mCFP/mYFP/KillerRed Antibody, Magnetic Beads

Product data:

Product Type:	Tag Antibodies
Product Name:	Mouse Monoclonal mGFP/mCFP/mYFP/KillerRed Antibody, Magnetic Beads
Clone Name:	OT12F6 (formerly 2F6)
Applications:	IP, Protein Purification
Recommended Dilution:	IP: Usage Dependent
Reactivity:	mGFP, mCFP, mYFP, KillerRed
Host:	Mouse
Clonality:	Monoclonal
Immunogen:	This antibody is raised against a recombinant protein of mGFP expressed in HEK293 cells. It recognizes mGFP, mCFP, mYFP and KillerRed proteins.
Isotype:	IgG1
Formulation:	PBS (pH7.4) containing 0.1% BSA, 0.02% NaN ₃ and 50% glycerol with final concentration 5mg/mL beads
Conjugation:	Magnetic Beads
Storage:	Shipped at -20°C or with ice packs, Upon delivery store at -20°C. Dilute in PBS(pH7.3) if necessary. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws.
Stability:	Stable for 1 year at -20°C from date of shipment.
Predicted Protein Size:	27.3 kDa



[View online »](#)