

Product datasheet for **TA183040S**

Mouse Monoclonal E2-Crimson Antibody, Magnetic Beads

Product data:

| | |
|-------------------------|--|
| Product Type: | Tag Antibodies |
| Product Name: | Mouse Monoclonal E2-Crimson Antibody, Magnetic Beads |
| Clone Name: | OTI9B3 (formerly 9B3) |
| Applications: | IP, Protein Purification |
| Recommended Dilution: | IP: Usage Dependent |
| Reactivity: | E2-Crimson |
| Host: | Mouse |
| Clonality: | Monoclonal |
| Immunogen: | This antibody is raised against a recombinant protein and recognizes the E2-Crimson protein. |
| Isotype: | IgG2a |
| Formulation: | PBS (pH7.4) containing 0.1% BSA, 0.02% NaN ₃ and 50% glycerol with final concentration 5mg/mL beads |
| Conjugation: | Magnetic Beads |
| Storage: | Shipped at -20°C or with ice packs, Upon delivery store at -20°C. Dilute in PBS(pH7.3) if necessary. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws. |
| Stability: | Stable for 1 year at -20°C from date of shipment. |
| Predicted Protein Size: | 26.2 kDa |



[View online »](#)