

Product datasheet for **TA183044**

Mouse Monoclonal Dendra2 Antibody, Magnetic Beads

Product data:

Product Type:	Tag Antibodies
Product Name:	Mouse Monoclonal Dendra2 Antibody, Magnetic Beads
Clone Name:	OT11G6 (formerly 1G6)
Applications:	IP, Protein Purification
Recommended Dilution:	IP: Usage Dependent
Reactivity:	Dendra2
Host:	Mouse
Clonality:	Monoclonal
Immunogen:	This antibody is raised against a recombinant protein of Dendra2 expressed in HEK293 cells.
Isotype:	IgG2b
Formulation:	PBS (pH7.4) containing 0.1% BSA, 0.02% NaN ₃ and 50% glycerol with final concentration 5mg/mL beads
Conjugation:	Magnetic Beads
Storage:	Shipped at -20°C or with ice packs, Upon delivery store at -20°C. Dilute in PBS(pH7.3) if necessary. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws.
Stability:	Stable for 1 year at -20°C from date of shipment.
Predicted Protein Size:	27.5 kDa



[View online »](#)