

## Product datasheet for **TA183048**

### Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody, DyLight 594

#### Product data:

Product Type:	Tag Antibodies
Product Name:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody, DyLight 594
Clone Name:	OTI4C5 (formerly 4C5)
Applications:	FC, IF
Recommended Dilution:	IF: 1:100, FC: 1:100
Reactivity:	DDK (FLAG tag): DYKDDDDK
Host:	Mouse
Clonality:	Monoclonal
Immunogen:	Anti-DDK monoclonal antibody is produced by immunizing mice with a synthetic peptide coupled to KLH
Isotype:	IgG2a
Formulation:	PBS (pH7.4) containing 0.1% BSA, 50% Glycerol with final concentration 500ug/mL
Concentration:	0.5mg/mL
Conjugation:	DyLight 594
Storage:	Store at -20°C. For maximum recovery of product, centrifuge the original vial after thawing and prior to opening the cap. Aliquot to avoid repeated freezing and thawing.
Stability:	Stable for 1 year at -20°C from date of shipment.
Synonyms:	FLAG TAG
Note:	This antibody doesn't recognize N- terminus DDK tag. Please use TA180144 for the detection of DDK tag at any location in Western blotting.



[View online »](#)