

Product datasheet for **TA183051**

Mouse monoclonal eGFP antibody, clone OT15A2, HRP conjugated

Product data:

| | |
|-----------------------|---|
| Product Type: | Tag Antibodies |
| Product Name: | Mouse monoclonal eGFP antibody, clone OT15A2, HRP conjugated |
| Clone Name: | OT15A2 (formerly 5A2) |
| Applications: | WB |
| Recommended Dilution: | WB: 1:1000 |
| Target: | eGFP |
| Host: | Mouse |
| Clonality: | Monoclonal |
| Immunogen: | A full length eGFP tagged recombinant protein expressed in HEK293T cell |
| Isotype: | IgG1 |
| Formulation: | PBS (pH7.4) containing 0.1% BSA, 50% Glycerol with final concentration 500ug/mL |
| Concentration: | 0.5mg/mL |
| Conjugation: | HRP |
| Storage: | Store at -20°C. For maximum recovery of product, centrifuge the original vial after thawing and prior to opening the cap. Aliquot to avoid repeated freezing and thawing. |
| Stability: | Stable for 1 year at -20°C from date of shipment. |



[View online »](#)