

Product datasheet for **TA183062**

HA mouse monoclonal antibody, clone CB051, Biotinylated

Product data:

Product Type:	Tag Antibodies
Product Name:	HA mouse monoclonal antibody, clone CB051, Biotinylated
Clone Name:	CB051
Applications:	FC, IF, WB
Recommended Dilution:	WB: 1:1000, IF: 1:100, FC: 1:100
Reactivity:	HA tag
Host:	Mouse
Clonality:	Monoclonal
Immunogen:	Monoclonal antibody is produced by immunizing mice with a synthetic peptide from influenza hemagglutinin coupled to KLH
Isotype:	IgG1
Formulation:	PBS (pH7.4) containing 0.1% BSA, 50% Glycerol with final concentration 500ug/mL
Concentration:	0.5mg/mL
Conjugation:	Biotin
Storage:	Store at -20°C. For maximum recovery of product, centrifuge the original vial after thawing and prior to opening the cap. Aliquot to avoid repeated freezing and thawing.
Stability:	Stable for 1 year at -20°C from date of shipment.



[View online »](#)