

Product datasheet for **TA183070**

HA Mouse Monoclonal Antibody, Clone CB051, Magnetic Beads

Product data:

Product Type:	Tag Antibodies
Product Name:	HA Mouse Monoclonal Antibody, Clone CB051, Magnetic Beads
Clone Name:	CB051
Applications:	IP, Protein Purification
Recommended Dilution:	IP, See Datasheet
Reactivity:	HA tag
Host:	Mouse
Clonality:	Monoclonal
Immunogen:	Monoclonal antibody is produced by immunizing mice with a synthetic peptide from influenza hemagglutinin coupled to KLH
Isotype:	IgG1
Formulation:	PBS (pH7.4) containing 0.1% BSA, 0.02% NaN ₃ and 50% glycerol with final concentration 5mg/mL beads
Conjugation:	Magnetic Beads
Storage:	Shipped at -20°C or with ice packs, Upon delivery store at -20°C. Dilute in PBS(pH7.3) if necessary. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws.
Stability:	Stable for 1 year at -20°C from date of shipment.



[View online »](#)