

## Product datasheet for **TA183081**

### Mouse Monoclonal tGFP Antibody, Clone OTI2H8, Agarose Beads

#### Product data:

Product Type:	Tag Antibodies
Product Name:	Mouse Monoclonal tGFP Antibody, Clone OTI2H8, Agarose Beads
Clone Name:	OTI2H8 (formerly 2H8)
Applications:	IP, Protein Purification
Recommended Dilution:	IP: See Datasheet
Target:	tGFP
Host:	Mouse
Clonality:	Monoclonal
Immunogen:	Anti-tGFP monoclonal antibody was produced by immunizing mice with purified tGFP protein expressed in HEK293T cell.
Isotype:	IgG2b
Formulation:	PBS (pH7.4) containing 0.1% BSA, 0.02% NaN <sub>3</sub> and 50% glycerol
Conjugation:	Agarose
Storage:	Shipped at -20°C or with ice packs, Upon delivery store at -20°C. Dilute in PBS(pH7.3) if necessary. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws.
Stability:	Stable for 1 year at -20°C from date of shipment.



[View online »](#)