

## **Product datasheet for TA190119**

## **BrdU Rat Monoclonal Antibody [Clone ID: OTI4A6]**

## **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI4A6

Applications: IHC

Recommended Dilution: IHC 1:150

**Reactivity:** Human

**Host:** Rat

**Isotype:** lgG2b

Clonality: Monoclonal

**Immunogen:** Bromodeoxyuridine coupled to keyhole limpet hemocyanin (KLH).

**Formulation:** PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

**Storage:** Store at -20°C as received.

**Stability:** Stable for 12 months from date of receipt.


**OriGene Technologies, Inc.** 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn


Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



## **Product images:**




Immunohistochemical staining of paraffinembedded colon tissue from BrdU injected mouse using anti-BrdU rat TA190119}



Immunohistochemical staining of paraffinembedded colon tissue from IDU injected mouse using anti-BrdU rat monoclonal antibody (TA190119).





Immunohistochemical staining of paraffinembedded colon tissue from CIDU injected mouse using anti-BrdU rat monoclonal antibody (TA190119).