

Product datasheet for TA190195

BrdU Rat Monoclonal Antibody [Clone ID: OTI4A3G8]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI4A3G8

Applications: IHC

Recommended Dilution: IHC 1:150

Reactivity: Human

Host: Rat

Isotype: lgG2b

Clonality: Monoclonal

Immunogen: Bromodeoxyuridine coupled to keyhole limpet hemocyanin (KLH).

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.


OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn


Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



Product images:




Immunohistochemical staining of paraffinembedded colon tissue from BrdU injected mouse using anti-BrdU rat monoclonal antibody (TA190195).



Immunohistochemical staining of paraffinembedded colon tissue from IDU injected mouse using anti-BrdU rat monoclonal antibody (TA190195).





Immunohistochemical staining of paraffinembedded colon tissue from CIDU injected mouse using anti-BrdU rat monoclonal antibody (TA190195).