

## **Product datasheet for TA301482**

## **OriGene Technologies, Inc.** 9620 Medical Center Drive, Ste 200

Rockville, MD 20850, US
Phone: +1-888-267-4436
https://www.origene.com
techsupport@origene.com
EU: info-de@origene.com
CN: techsupport@origene.cn

## **BMI1 Rabbit Polyclonal Antibody**

**Product data:** 

**Product Type:** Primary Antibodies

**Applications:** ICC/IF, IHC, WB

**Recommended Dilution:** Immunohistochemistry, Western Blot: 2 ug/ml, Immunocytochemistry/ Immunofluorescence:

1:50-1:200

Reactivity: Human, Mouse

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

Immunogen: A synthetic peptide made to an internal region (within residues 1-100) of the human Bmi1

protein. [Swiss-Prot# P35226]

Formulation: Tris-glycine, 150mM NaCl and 0.05% sodium azide

**Concentration:** lot specific

**Purification:** Immunogen affinity purified

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Gene Name: BMI1 proto-oncogene, polycomb ring finger

Database Link: NP 005171

Entrez Gene 12151 MouseEntrez Gene 648 Human

P35226

**Background:** The Bmi1 oncogene has been shown to induce telomerase activity and immortalize human

mammary epithelial cells. Bmi1 extends the replicative life span of human fibroblasts by suppressing the p16-dependent senescence pathway. Bmi1 is also important in maintaining

stem cells' ability to self renew. It has also been termed an oncogene with aberrant expression in various carcinomas. It has been shown to be a negative regulator of p19Arf,

Ink4a (p14Arf) and tp53 known tumor repressors. Ablation of Bmi1 is required for short term

survival of cancer cells.

Synonyms: BMI1; FLVI2; PCGF4; RNF51







**Protein Families:** 

Druggable Genome, ES Cell Differentiation/IPS, Transcription Factors

## **Product images:**



Immunocytochemistry/Immunofluorescence: Bmi1 Antibody TA301482 - BMI1 antibody was tested in HeLa cells with Dylight 488 (green). Nuclei and alpha-tubulin were counterstained with DAPI (blue) and Dylight 550 (red).



Western Blot: Bmi1 Antibody TA301482 -Detection of Bmi1 protein in skeletal muscle