

Product datasheet for TA302127

CaMKI gamma (CAMK1G) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 1:1000, IHC: 1:50~100
Reactivity:	Human
Host:	Rabbit
Isotype:	Ig
Clonality:	Polyclonal
Immunogen:	This CAMK1G (CaMKI gamma) antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 420-450 amino acids from the C-terminal region of human CAMK1G (CaMKI gamma).
Formulation:	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Concentration:	lot specific
Purification:	This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	53128 Da
Gene Name:	calcium/calmodulin dependent protein kinase IG
Database Link:	NP_065172 Entrez Gene 57172 Human
Background:	Ca ²⁺ /calmodulin-dependent protein kinase I (CaMKI) constitutes a family of closely related isoforms (alpha, beta and gamma). CLICK-III/CaMKIgamma is a novel membrane-anchored neuronal Ca ²⁺ /calmodulin-dependent protein kinase. AMKlgamma is abundant in neurons, particularly in the amygdala and ventromedial hypothalamus. Like the other CaMKI isoforms, full activation of CLICK-III/CaMKIgamma requires both Ca(2+)/CaM and phosphorylation by CaMKK.
Synonyms:	CLICK3; CLICKIII; dj272L16.1; VWS1
Protein Families:	Druggable Genome, Protein Kinase




[View online »](#)


Product images:



Western blot analysis of anti-CAMK 1G Pab (Cat. #TA302127) in NCI-H460 cell lysate. CAMK 1G (arrow) was detected using purified Pab. Secondary HRP-anti-rabbit was used for signal visualization with chemiluminescence.



Western blot analysis of CAMK1G (arrow) using rabbit polyclonal CAMK1G Antibody (C-term) ([RB01249]). 293 cell lysates (2 ug/lane) either nontransfected (Lane 1) or transiently transfected with the CAMK1G gene (Lane 2) (Origene Technologies).



Formalin-fixed and paraffin-embedded human cancer tissue reacted with the primary antibody, which was peroxidase-conjugated to the secondary antibody, followed by AEC staining. This data demonstrates the use of this antibody for immunohistochemistry; clinical relevance has not been evaluated. BC = breast carcinoma; HC = hepatocarcinoma.