

## Product datasheet for TA302551

### Nogo Receptor (RTN4R) Goat Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	ELISA: 1:16,000. WB: 0.3-1 µg/ml.
Reactivity:	Human (Expected from sequence similarity: Mouse, Rat, Dog, Pig)
Host:	Goat
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Peptide with sequence C-LQALPDDTFRDLGN, from the internal region of the protein sequence according to NP_075380.1.
Formulation:	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Concentration:	lot specific
Purification:	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide. Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	54335 Da
Gene Name:	reticulon 4 receptor
Database Link:	<a href="#">NP_075380 Entrez Gene 65079 Mouse</a> <a href="#">Entrez Gene 113912 Rat</a> <a href="#">Entrez Gene 486414 Dog</a> <a href="#">Entrez Gene 65078 Human</a>
Background:	This gene encodes the receptor for reticulon 4, oligodendrocyte myelin glycoprotein and myelin-associated glycoprotein. This receptor mediates axonal growth inhibition and may play a role in regulating axonal regeneration and plasticity in the adult central nervous system. [provided by RefSeq]
Synonyms:	NGR; NOGOR
Protein Families:	Druggable Genome



[View online »](#)

## Product images:



TA302551 (0.3ug/ml) staining of Human Brain lysate (35ug protein in RIPA buffer) with (B) and without (A) blocking with the immunising peptide. Primary incubation was 1 hour. Detected by chemiluminescence.