

Product datasheet for **TA303552**

IRS1 Rabbit Monoclonal Antibody [Clone ID: EP260Y]

Product data:

Product Type:	Primary Antibodies
Clone Name:	EP260Y
Recommended Dilution:	WB: 1:500000
Reactivity:	Mouse, Rat, Human
Modifications:	Phospho-specific
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	A synthetic peptide corresponding to residues near the C-terminus of human IRS-1 was used as immunogen.
Formulation:	PBS 49%, Sodium azide 0.01%, Glycerol 50%, BSA 0.05%
Purification:	Tissue culture supernatant
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	132 kDa
Gene Name:	insulin receptor substrate 1
Database Link:	NP_005535 Entrez Gene 16367 Mouse Entrez Gene 25467 Rat Entrez Gene 3667 Human P35568
Synonyms:	HIRS-1
Note:	Is unsuitable for Flow Cyt, ICC/IF, IHC-P or IP.
Protein Families:	Druggable Genome
Protein Pathways:	Adipocytokine signaling pathway, Insulin signaling pathway, Neurotrophin signaling pathway, Type II diabetes mellitus



[View online »](#)