

Product datasheet for TA305818

http tech

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US

OriGene Technologies, Inc.

Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

MAML1 Goat Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: WB

Recommended Dilution: WB: 0.1-0.3ug/ml.

Reactivity: Human (Expected from sequence similarity: Mouse)

Host: Goat Isotype: IgG

Clonality: Polyclonal

Immunogen: Peptide with sequence VLPTCPMAEFALPR, from the N Terminus of the protein sequence

according to NP_055572.1.

Formulation: 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin

Concentration: lot specific

Purification: Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity

chromatography using the immunizing peptide. Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. Aliquot and store at -20C. Minimize

freezing and thawing.

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Gene Name: mastermind like transcriptional coactivator 1

Database Link: NP 055572

Entrez Gene 103806 MouseEntrez Gene 9794 Human

Q92585

Background: This protein is the human homolog of mastermind, a Drosophila protein that plays a role in

the Notch signaling pathway involved in cell-fate determination. There is in vitro evidence that the human homolog forms a complex with the intracellular portion of human Notch receptors and can increase expression of a Notch-induced gene. This evidence supports its proposed function as a transcriptional co-activator in the Notch signaling pathway. [provided

by RefSeq1





MAML1 Goat Polyclonal Antibody - TA305818

Synonyms: Mam-1; Mam1

Protein Pathways: Notch signaling pathway

Product images:

250kDa 150kDa 100kDa 75kDa 50kDa 37kDa

37kDa 25kDa 20kDa 15kDa TA305818 (0.1ug/ml) staining of Human Skeletal Muscle lysate (35ug protein in RIPA buffer). Primary incubation was 1 hour. Detected by

chemiluminescence.