

Product datasheet for **TA305985**

Apoptosis repressor with CARD (NOL3) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	WB: 1 - 2 ug/mL, IHC: 5 ug/mL
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	ARC antibody was raised against a peptide corresponding to amino acids 191 to 208 of human origin.
Formulation:	PBS containing 0.02% sodium azide.
Concentration:	1ug/ul
Purification:	Affinity chromatography purified via peptide column
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	nucleolar protein 3
Database Link:	NP_003937 Entrez Gene 8996 Human O60936
Background:	Apoptosis is regulated by death domain (DD) and/or caspase recruitment domain (CARD) containing molecules and a caspase family of proteases. CARD domain containing cell death regulators include RAIDD, Apaf-1, caspase-9, and caspase-2. A novel CARD domain containing protein was recently identified and designated ARC for apoptosis repressor with CARD (1). ARC interacts with caspase-2 and -8 and inhibits enzymatic activity of caspase-8. ARC suppresses apoptosis induced by cell death adapters FADD and TRADD and by cell death receptors Fas, TNFR-1 and DR3. The messenger RNA of ARC is primarily expressed in skeletal muscle and cardiac tissue (1).
Synonyms:	ARC; FCM; MYP; NOP; NOP30




[View online »](#)


Product images:



Western blot analysis of ARC in HeLa (H), KB (K), and A549 (A) whole cell lysates with ARC antibody at 1:1000 dilution. 1st Below: Immunohistochemistry of ARC in human heart tissue at 5 ug/ml. 2nd Below: Immunohistochemistry of ARC in human skeletal muscle with ARC antibody at 5 ug/ml.



Immunohistochemistry of ARC in human skeletal muscle with ARC antibody at 5 ug/ml.



Immunohistochemistry of ARC in human heart tissue with ARC antibody at 5 ug/mL.