

## **Product datasheet for TA306345**

## PAK2 Rabbit Polyclonal Antibody

## **Product data:**

**Product Type:** Primary Antibodies

**Applications:** IF, IHC, WB

Recommended Dilution: WB: 0.5 - 1 ug/mL, ICC: 10 ug/mL, IF: 20 ug/mL

Reactivity: Human, Mouse, Rat

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

**Immunogen:** PAK2 antibody was raised against a 14 amino acid peptide from near the carboxy terminus of

human PAK2.

**Formulation:** PBS containing 0.02% sodium azide.

Concentration: 1ug/ul

**Purification:** Affinity chromatography purified via peptide column

**Conjugation:** Unconjugated

**Storage:** Store at -20°C as received.

**Stability:** Stable for 12 months from date of receipt.

**Gene Name:** p21 (RAC1) activated kinase 2

Database Link: NP 002568

Entrez Gene 29432 RatEntrez Gene 224105 MouseEntrez Gene 5062 Human

Q13177



**OriGene Technologies, Inc.** 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



Background:

The p21-activated kinases (PAKs) are serine-threonine kinases that bind to the active forms of Cdc42 and Rac. They are divided into two groups, the first of which include PAK1, 2 and 3, and can be activated by Cdc42/Rac binding. Group 1 PAKs contain an autoinhibitory domain whose activity is regulated by Cdc42/Rac binding. The group 1 PAKs are known to be involved in cellular processes such as gene transcription, apoptosis, and cell morphology and motility. Much less is known about the second group, which includes PAK4, 5 and 6, and are not activated by Cdc42/Rac binding. Of the six PAK proteins, only PAK2 is ubiquitously expressed and cleaved by caspase-3. This cleavage removes the amino-terminal regulatory domain and generates a constitutively active kinase fragment. Recent experiments have shown that following cleavage, the active fragment is myristoylated and directed to the plasma membrane and membrane ruffles where it promotes cell death via increased signaling through the c-Jun N-terminal kinase pathway, but without compromising mitochondrial integrity.


Synonyms: PAK65; PAKgamma

**Protein Families:** Druggable Genome, Protein Kinase


**Protein Pathways:** Axon guidance, ErbB signaling pathway, Focal adhesion, MAPK signaling pathway, Regulation

of actin cytoskeleton, Renal cell carcinoma, T cell receptor signaling pathway

## **Product images:**




Western blot analysis of PAK2 in Jurkat lysate with PAK2 antibody at (A) 0.5 and (B) 1 ug/ml.



Immunohistochemistry of PAK2 in mouse spleen tissue with PAK2 antibody at 10 ug/ml.





Immunofluorescence of PAK2 in Mouse Spleen cells with PAK2 antibody at 20 ug/mL.