

Product datasheet for **TA306490**

FNBP1L Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC, WB
Recommended Dilution:	WB: 1 ug/mL, ICC: 2.5 ug/mL, IF: 20 ug/mL
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	TOCA-1 antibody was raised against a 17 amino acid peptide near the carboxy terminus of the human TOCA-1.
Formulation:	PBS containing 0.02% sodium azide.
Concentration:	1ug/ul
Purification:	Affinity chromatography purified via peptide column
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	formin binding protein 1 like
Database Link:	NP_001020119 Entrez Gene 214459 Mouse Entrez Gene 310839 Rat Entrez Gene 54874 Human Q5TON5




[View online »](#)

Background:


Actin reorganization is important for the regulation of neuronal morphology. A protein involved in this process, the transducer of cdc42-dependent actin assembly 1 (TOCA-1) protein, a member of the evolutionarily conserved pombe CDC15 homology (PCH) protein family, is an essential component of the Cdc42 pathway. TOCA-1 binds both N-WASP and Cdc42 and is essential for Cdc42- and PIP2-induced actin polymerization, suggesting that TOCA-1 mediates Cdc42-dependent actin nucleation by activating the N-WASP-WIP complex. Decreased expression of TOCA-1 significantly enhances neurite elongation in PC-12 cells; its overexpression in the same cells suppresses neurite elongation. It has been suggested that TOCA-1 negatively regulates axon branching by regulating membrane trafficking by regulating membrane trafficking through the F-BAR/EFC domain. Multiple isoforms of TOCA-1 are known to exist.

Synonyms:


C1orf39; TOCA1

Product images:

Western blot analysis of TOCA-1 in mouse brain tissue lysate with in (A) the absence and (B) the presence of blocking peptide with TOCA-1 antibody at 1 ug/ml.



Immunohistochemistry of TOCA-1 in human brain tissue with TOCA-1 antibody at 2.5 ug/ml.



Immunofluorescence of TOCA-1 in Human Brain cells with TOCA-1 antibody at 20 ug/mL.