

Product datasheet for **TA308255**

Syntaxin 4 (STX4) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB:1:500-1:3000
Reactivity:	Human (Predicted: Rat, Bovine)
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant fragment corresponding to a region within amino acids 92 and 297 of Syntaxin 4 (Uniprot ID#Q12846)
Formulation:	0.1M Tris, 0.1M Glycine, 20% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.
Concentration:	lot specific
Purification:	Purified by antigen-affinity chromatography.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	34 kDa
Gene Name:	syntaxin 4
Database Link:	NP_004595 Entrez Gene 81803 Rat Entrez Gene 6810 Human Q12846
Synonyms:	p35-2; STX4A
Note:	Seq homology of immunogen across species: Rat (97%), Bovine (95%)
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	SNARE interactions in vesicular transport



[View online »](#)

Product images:



Sample (30 ug of whole cell lysate). A: A431. 10% SDS PAGE. TA308255 diluted at 1:3000.