

Product datasheet for TA308320

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com

OriGene Technologies, Inc.

EU: info-de@origene.com CN: techsupport@origene.cn

Ataxin 3 (ATXN3) Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: IHC, WB

Recommended Dilution: IHC:1:100-1:1000; WB:1:500-1:10000

Reactivity: Human, Mouse (Predicted: Pig, Rat, Bovine)

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

Immunogen: Recombinant fragment corresponding to a region within amino acids 1 and 245 of Ataxin 3

(Uniprot ID#P54252)

Formulation: 1XPBS, 1% BSA, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative.

Concentration: lot specific

Purification: Purified by antigen-affinity chromatography.

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 42 kDa **Gene Name:** ataxin 3

Database Link: NP 001121169

Entrez Gene 60331 RatEntrez Gene 110616 MouseEntrez Gene 4287 Human

P54252

Background: Machado-Joseph disease, also known as spinocerebellar ataxia-3, is an autosomal dominant

neurologic disorder. The protein encoded by this gene contains (CAG)n repeats in the coding region, and the expansion of these repeats from the normal 13-36 to 68-79 is the cause of Machado-Joseph disease. There is a negative correlation between the age of onset and CAG repeat numbers. Alternatively spliced transcript variants encoding different isoforms have

been described for this gene. [provided by RefSeq]

Synonyms: AT3; ATX3; JOS; MJD; MJD1; SCA3







Note: Seq homology of immunogen across species: Pig (99%), Rat (96%), Bovine (98%)

Protein Families: Druggable Genome, Transcription Factors

Product images:




Ataxin 3 antibody detects Ataxin 3 protein by western blot analysis. Mouse tissue extracts (50 ug) was separated by 10 % SDS-PAGE, and blotted with Ataxin 3 antibody (TA308320) diluted by 1:500



Ataxin 3 antibody detects Ataxin 3 protein by western blot analysis. 293T mock transfected (-) or transfected with a construct expressing human Ataxin 3 (+). 10 % SDS-PAGE. Ataxin 3 antibody (TA308320) dilution: 1:5000





Immunohistochemical analysis of paraffinembedded H661 xenograft, using Ataxin 3 (TA308320) antibody at 1:500 dilution.