

## Product datasheet for **TA308329**

### **CENTB1 (ACAP1) Rabbit Polyclonal Antibody**

#### **Product data:**

<b>Product Type:</b>	Primary Antibodies
<b>Reactivity:</b>	Human (Predicted: Mouse, Rat, Bovine)
<b>Host:</b>	Rabbit
<b>Isotype:</b>	IgG
<b>Clonality:</b>	Polyclonal
<b>Immunogen:</b>	Synthetic peptide corresponding to a region within amino acids 678 and 740 of Centaurin beta 1 (Uniprot ID#Q15027)
<b>Formulation:</b>	0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.
<b>Concentration:</b>	lot specific
<b>Purification:</b>	Purified by antigen-affinity chromatography.
<b>Conjugation:</b>	Unconjugated
<b>Storage:</b>	Store at -20°C as received.
<b>Stability:</b>	Stable for 12 months from date of receipt.
<b>Predicted Protein Size:</b>	82 kDa
<b>Gene Name:</b>	ArfGAP with coiled-coil, ankyrin repeat and PH domains 1
<b>Database Link:</b>	<a href="#">NP_055531</a> <a href="#">Entrez Gene 216859 MouseEntrez Gene 287443 RatEntrez Gene 9744 Human Q15027</a>
<b>Synonyms:</b>	CENTB1
<b>Note:</b>	Seq homology of immunogen across species: Mouse (100%), Rat (100%), Bovine (100%)
<b>Protein Pathways:</b>	Endocytosis



[View online »](#)