

Product datasheet for TA308348

GALR2 Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: IHC, WB

Recommended Dilution: IHC:1:100-1:1000; WB:1:500-1:3000

Reactivity: Human (Predicted: Bovine, Chicken, Mouse)

Host: Rabbit Isotype: lgG

Polyclonal Clonality:

Immunogen: Synthetic peptide corresponding to a region within amino acids 243 and 336 of Galanin

Receptor 2 (Uniprot ID#O43603)

Formulation: 1XPBS, 50% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.

Concentration: lot specific

Purification: Purified by antigen-affinity chromatography.

Conjugation: Unconjugated

Store at -20°C as received. Storage:

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 42 kDa

Gene Name: galanin receptor 2

Database Link: NP 003848

Entrez Gene 14428 MouseEntrez Gene 8811 Human

O43603



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



Background:


Galanin is an important neuromodulator present in the brain, gastrointestinal system, and hypothalamopituitary axis. It is a 30-amino acid non-C-terminally amidated peptide that potently stimulates growth hormone secretion, inhibits cardiac vagal slowing of heart rate, abolishes sinus arrhythmia, and inhibits postprandial gastrointestinal motility. The actions of galanin are mediated through interaction with specific membrane receptors that are members of the 7-transmembrane family of G protein-coupled receptors. GALR2 interacts with the N-terminal residues of the galanin peptide. The primary signaling mechanism for GALR2 is through the phospholipase C/protein kinase C pathway (via Gq), in contrast to GALR1, which communicates its intracellular signal by inhibition of adenylyl cyclase through Gi. However, it has been demonstrated that GALR2 couples efficiently to both the Gq and Gi proteins to simultaneously activate 2 independent signal transduction pathways. [provided by RefSeq]

Synonyms: GAL2-R; GALNR2; GALR-2

Note: Seq homology of immunogen across species: Mouse (92%), Chicken (86%), Bovine (93%)


Protein Families: Druggable Genome, Transmembrane
Protein Pathways: Neuroactive ligand-receptor interaction

Product images:




Sample (30 ug of whole cell lysate). A:A431. 12% SDS PAGE. TA308348 diluted at 1:500





Galanin Receptor 2 antibody detects GALR2 protein at membrane on Cal27 xenograft by immunohistochemical analysis. Sample: Paraffinembedded Cal27 xenograft. Galanin Receptor 2 antibody (TA308348) dilution: 1:500.



Galanin Receptor 2 antibody detects Galanin Receptor 2 protein at membrane and cytosol on rat fore brain by immunohistochemical analysis. Sample: Paraffin-embedded rat fore brain. Galanin Receptor 2 antibody (TA308348) dilution: 1:500