

# Product datasheet for TA308364

# CD98 (SLC3A2) Rabbit Polyclonal Antibody

## **Product data:**

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	IHC:1:100-1:1000; WB:1:500-1:3000
Reactivity:	Human (Predicted: Rabbit, Cow)
Host:	Rabbit
lsotype:	lgG
Clonality:	Polyclonal
Immunogen:	Recombinant fragment corresponding to a region within amino acids 89 and 433 of CD98 (Uniprot ID#P08195)
Formulation:	1XPBS, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative.
Concentration:	lot specific
Purification:	Purified by antigen-affinity chromatography.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	71 kDa
Gene Name:	solute carrier family 3 member 2
Database Link:	<u>NP_002385</u> <u>Entrez Gene 6520 Human</u> <u>P08195</u>
Background:	This gene is a member of the solute carrier family and encodes a cell surface, transmembrane protein with an alpha amylase domain. The protein exists as the heavy chain of a heterodimer, covalently bound through di-sulfide bonds to one of several possible light chains. It associates with integrins and mediates integrin-dependent signaling related to normal cell growth and tumorigenesis. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq]
Synonyms:	4F2; 4F2HC; 4T2HC; CD98; CD98HC; MDU1; NACAE



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

## **CD98** (SLC3A2) Rabbit Polyclonal Antibody – TA308364


#### Seq homology of immunogen across species: Rabbit (83%), cow (81%)

**Protein Families:** 


Note:

: Transmembrane

# **Product images:**



Sample (30 ug of whole cell lysate). A: Molt-4. 7.5% SDS PAGE. TA308364 diluted at 1:1000.



Immunohistochemical analysis of paraffinembedded SCC15 xenograft, using CD98 (TA308364) antibody at 1:100 dilution.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US