

Product datasheet for TA308405

IPMK Rabbit Polyclonal Antibody

Product data:

OriGene Technologies, Inc.


9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product Type:	Primary Antibodies
Applications:	IF, WB
Recommended Dilution:	ICC/IF:1:100-1:1000; WB:1:500-1:3000
Reactivity:	Human, Mouse (Predicted: Rat, Bovine)
Host:	Rabbit
lsotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant fragment corresponding to a region within amino acids 26 and 295 of IPMK (Uniprot ID#Q8NFU5)
Formulation:	0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.
Purification:	Purified by antigen-affinity chromatography.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	47 kDa
Gene Name:	inositol polyphosphate multikinase
Database Link:	<u>NP_689416</u> Entrez Gene 69718 MouseEntrez Gene 171458 RatEntrez Gene 253430 Human <u>Q8NFU5</u>
Synonyms:	Impk; inositol polyphosphate kinase 2; inositol polyphosphate multikinase; Ipk2
Note:	Seq homology of immunogen across species: Rat (90%), Bovine (92%)
Protein Pathways:	Inositol phosphate metabolism




This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

Product images:




Sample (50 ug of whole cell lysate). A: Mouse brain. 10% SDS PAGE. TA308405 diluted at 1:1000.



Sample (30 ug of whole cell lysate). A: H1299. 10% SDS PAGE. IPMK antibody. TA308405 diluted at 1:1000.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US





Immunofluorescence analysis of paraformaldehyde-fixed HeLa, using IPMK (TA308405) antibody at 1:200 dilution.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US