

## Product datasheet for **TA308760**

### Cullin 5 (CUL5) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Reactivity:	Human (Predicted: Mouse, Feline, Rabbit)
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Recombinant fragment corresponding to a region within amino acids 113 and 583 of Cullin 5 (Uniprot ID#Q93034)
Formulation:	0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.
Concentration:	lot specific
Purification:	Purified by antigen-affinity chromatography.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	91 kDa
Gene Name:	cullin 5
Database Link:	<a href="#">NP_003469</a> <a href="#">Entrez Gene 75717 Mouse</a> <a href="#">Entrez Gene 8065 Human</a> <a href="#">Q93034</a>
Synonyms:	VACM-1; VACM1
Note:	Seq homology of immunogen across species: Human (100%), Mouse (100%), Cat (100%), Rabbit (100%)
Protein Families:	Druggable Genome
Protein Pathways:	Ubiquitin mediated proteolysis



[View online »](#)