

Product datasheet for TA308837

Actin (ACTA1) Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: IHC, WB

Recommended Dilution: WB 1:1000 - 1:10000, IHC:1:100 - 1:1000

Reactivity: Human, Mouse (Predicted: Rat, Pig, Chicken, Bovine)

Host: Rabbit Isotype: lgG

Clonality: Polyclonal

Recombinant protein fragment contain a sequence corresponding to a region within amino Immunogen:

acids 1 and 377 of alpha Actin (skeletal muscle)

Formulation: 1XPBS, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative.

Concentration: lot specific

Purification: Purified by antigen-affinity chromatography.

Conjugation: Unconjugated

Store at -20°C as received. Storage:

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 42 kDa

Gene Name: actin, alpha 1, skeletal muscle

Database Link: NP 001091

Entrez Gene 11459 MouseEntrez Gene 29437 RatEntrez Gene 58 Human

P68133

Background: The product encoded by this gene belongs to the actin family of proteins, which are highly

> conserved proteins that play a role in cell motility, structure and integrity. Alpha, beta and gamma actin isoforms have been identified, with alpha actins being a major constituent of the contractile apparatus, while beta and gamma actins are involved in the regulation of cell motility. This actin is an alpha actin that is found in skeletal muscle. Mutations in this gene cause nemaline myopathy type 3, congenital myopathy with excess of thin myofilaments, congenital myopathy with cores, and congenital myopathy with fiber-type disproportion,

diseases that lead to muscle fiber defects. [provided by RefSeq]



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com




Synonyms: ACTA; ASMA; CFTD; CFTD1; CFTDM; MPFD; NEM1; NEM2; NEM3; SHPM

Note: Seq homology of immunogen across species: Rat (100%), Pig (100%), Chicken (100%), Bovine


(100%

Protein Families: Stem cell - Pluripotency

Product images:




Sample (30 ug of whole cell lysate). A: NIH-3T3. B: JC. C: BCL-1. D: C2C12. 10% SDS PAGE. TA308837 diluted at 1:5000.



Sample (30 ug of whole cell lysate). A:Raji. 10% SDS PAGE. TA308837 diluted at 1:1000





Immunohistochemical analysis of paraffinembedded Skeletal Muscle, using alpha Actin (skeletal muscle) (TA308837) antibody (10 ug/ml).