

## Product datasheet for **TA308911**

### ISLR Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.
Concentration:	lot specific
Purification:	Purified by antigen-affinity chromatography.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	46 kDa
Gene Name:	immunoglobulin superfamily containing leucine-rich repeat
Database Link:	<a href="#">NP_958934</a> <a href="#">Entrez Gene 3671 Human</a> <a href="#">O14498</a>
Synonyms:	HsT17563



[View online »](#)