

## Product datasheet for **TA308987**

### beta III Tubulin (TUBB3) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	0.1M Tris, 0.1M Glycine, 10% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.
Concentration:	lot specific
Purification:	Purified by antigen-affinity chromatography.
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	50 kDa
Gene Name:	tubulin beta 3 class III
Database Link:	<a href="#">NP_006077</a> <a href="#">Entrez Gene 22152 Mouse</a> <a href="#">Entrez Gene 10381 Human</a> <a href="#">Q13509</a>
Synonyms:	beta-4; CDCBM; CDCBM1; CFEOM3; CFEOM3A; FEOM3; TUBB4
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS
Protein Pathways:	Gap junction, Pathogenic Escherichia coli infection



[View online »](#)