

Product datasheet for **TA309098**

CNPase (CNP) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	WB: 1:100, IHC: 1:1000
Reactivity:	Mouse, Human, Rabbit, Rat, Sheep
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Endogenous rabbit 2,3 cyclic nucleotide-3-phospho-diesterase
Formulation:	Neat serum
Purification:	Neat serum
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	46 kDa
Gene Name:	2',3'-cyclic nucleotide 3' phosphodiesterase
Database Link:	NP_149124 Entrez Gene 12799 Mouse Entrez Gene 25275 Rat Entrez Gene 1267 Human P09543
Synonyms:	CNP1



[View online »](#)

Product images:



Western blot of rat brain lysate showing the specific immunolabeling of the ~ 46k CNP protein.